

freeway

flow

Shower Cradle



User Manual

Contents

1	Safety Instructions and Warnings	3
1.1	Introduction.....	3
1.2	Manufacture.....	3
1.3	European Authorised Representative	3
1.4	Symbols Used.....	4
1.5	Contraindications / Limitations.....	4
1.6	Intended Use.....	4
1.7	Warning Notes	6
1.8	Operating Environment.....	6
2	Components/Key Parts.....	7
3	Accessories and Spare Parts.....	8
3.1	Accessories	8
3.2	Spare Parts.....	9
4	Assembly Instructions	10
4.1	Unpacking	10
5	Operating Instructions	11
5.1	Operating the Tilt Lever	11
5.2	Operating the Braked Castors.....	11
5.3	Attaching & Detaching the Shower Cradle Mesh.....	12
5.4	Manually Adjusting the Product Height.....	14
5.5	Using the Commode Pan.....	15
5.6	Positioning the Head Cushion.....	15
5.7	Accessories	16
6	Technical Specification.....	17
6.1	Shower Cradle Dimensions.....	17
6.2	Specifications.....	18
6.3	Expected Product Lifetime.....	18
7	Environmental - Storage and Operating Conditions	18
7.1	Normal Operating Conditions.....	18
7.2	Shipping and Storage Conditions.....	18
8	Disposal.....	19
9	Fault Finding	20
10	General Inspection and Maintenance.....	20
10.1	Service.....	20
10.2	Inspection.....	21
11	Cleaning.....	21
11.1	Washing & Drying of the Mesh Fabric.....	21
11.2	General Cleaning	21
11.3	Disinfecting (if necessary)	21
12	Warranty.....	23
13	Service Record History.....	24

1 Safety Instructions and Warnings

1.1 Introduction



Please read and understand this manual in its entirety before using your Freeway Flow Shower Cradle.

The information in this manual is important for the safety of anyone near the shower cradle and must be read and understood to help prevent injuries. It is also crucial to the proper operation and maintenance of the shower cradle.

Should any questions arise from reviewing this manual, contact your local authorised representative.

Failure to comply with warnings in this manual may result in; injury to the operator and/or client and/or damage to the shower cradle or related components.

If, during the use of this device or because of its use, a serious incident has occurred please report it to the manufacturer and to your national authority.

Store this manual with the documents included with the shower cradle. Contents of this manual are subject to change without prior written notice.



Do not attempt to use this equipment without first understanding the contents of this manual.



Unauthorised modifications on this product may affect its safety. The manufacturer will not be held responsible for any accident, incident or deficiencies of performance that occur because of any unauthorised modification to its products.

1.2 Manufacture

The shower cradle is manufactured at the address below:



Prism Medical UK

Unit 1, Tir Llwyd Industrial Estate, St Asaph Avenue, Kinmel Bay, Conwy, LL18 5JZ

Telephone number: 01924 840100

1.3 European Authorised Representative

The address of the European Authorised Representative for this product:



European Healthcare & Device Solutions (Ireland) Ltd

Stratton House, Bishopstown Road, Cork, Ireland, T12 Y9TC

Telephone number: +353(86)2280846

1.4 Symbols Used

The Table below includes all symbols from BS EN ISO 15223-1:2021 that can be found in this manual and on the Product and what they represent. Refer to this Table when you are unsure of what a symbol represents.











	Consult instructions before use		Caution – see instructions for use
	Manufacturer		Date of manufacture
	Catalogue number		EC Authorised Representative
	Serial number		Humidity range
	Temperature range		Packaging indicator – Keep dry

Table 1

1.5 Contraindications / Limitations

There are no known “contraindications” associated with the usage of the shower cradle, provided it is used as per the manufacturer’s recommendations and guidelines. However, it is recommended that a client specific assessment is completed by a trained and knowledgeable health care professional to determine the method of transfer and use.

The manufacturer does not recommend a required number of care givers for the use of our products. This information and recommendation can only be provided after a thorough personalised, case specific assessment, as there are many factors that can influence these decisions.

1.6 Intended Use

For internal use only.

The shower cradle, which is designed to provide safe, and comfortable seating and showering of a user. The manufacturer always recommends training before use of our products, for more details on the training that can be provided please contact your local authorised dealer.

All four castors come as standard with independent locking brakes which, makes it easy for a carer to manoeuvre, and provides maximum stability and safety for the user. The castors can also be manually height adjusted to suit carers of different heights.



The shower cradle has been tested at our manufacturing facility to a safe working load of 130 kg, this load must not be exceeded.
The shower cradle is used under instruction, and the operation of the aid is undertaken by a trained carer.

The shower cradle mesh is a specially designed fabric that attaches to the shower cradle by a system of buckles. The mesh can be removed for cleaning/washing and drying and is expected to be replaced during the lifetime of the metalwork. Please read and understand the section in this manual, for the correct attachment and detachment of the mesh to ensure the safety of the user is not compromised.

- If you are unsure on the correct use of this product, please contact the manufacturer or a professional for further information or training.
- The user must only be manoeuvred in the transport position, this is when the cradle is fully tilted up.
- Do not attempt to tilt the shower cradle when there is no user present, as the gas struts are only designed to operate with weight present.
- The shower cradle is intended must only be used on a flat and levelled surface.
- Under no circumstance should the shower cradle be put in control of a person who has not been properly trained in the use and care of this equipment. Failure to adhere to this warning may result in serious injury to the operator and/or the user.
- The shower cradle, and associated accessories are not toys. Do not use it for unsafe practices. Do not allow children to play with the shower cradle or any of its components. The shower cradle should not be used for any practice except its intended use.
- In facilities where more than one operator will be responsible for using the shower cradle, it is important that all such members be trained on the shower cradle prior to use. A training program should be established by the facility to acquaint new operators with this equipment.
- Your guarantee is void if persons unauthorised by the manufacturer perform work on the shower cradle.
- To maintain optimum function, the shower cradle should be inspected and maintained on a regular basis. See section 'General Inspection, Maintenance and Cleaning' within this user manual.
- A full list of approved accessories and spare parts are listed in this user manual.
- The shower cradle and the associated accessories are intended only for showering and transferring of a person. The manufacturer will not be responsible for any damage caused by the misuse, neglect, or purposeful destruction of the equipment and/or its associated components.
- Any accessories used with the shower cradle including the mesh fabric, should be checked before each use to ensure that they are in good working order. Check for signs of wear and ensure that all labels are legible. Report any unusual wear to your local authorised dealer.
- The shower cradle and its associated accessories are certified to a maximum load of 130 kg. Do not exceed the maximum rated load of any of the components.
- Ensure that a clear space is maintained around the shower cradle. Before using the shower cradle always check for and move away any obstacles.
- Never leave a user unattended in the shower cradle.
- The shower cradle is not intended to be used as a seating product for long periods of time.



You may need to seek specialist advice on how to assist some people with specific moving and handling needs. Sources of advice include, but is not limited to, professional bodies and organisations, occupational therapists, physiotherapists, manual handling advisers and ergonomists with experience in health and social care.

1.7 Warning Notes



The shower cradle has been designed, tested, and manufactured to a high standard. This does not mean that it can be used without care. All operators and carers should read the user manual and understand this warning section.

1. All shower cradles are less stable on sloping surfaces. A 5° slope is the maximum permitted.
2. All shower cradles are dangerous to the person being transferred when used recklessly or pushed at speed.
3. All shower cradles are less stable when the load is moving.
4. All shower cradles have less grip when pushed on a wet surface, where this product is intended for use.

1.8 Operating Environment

The shower cradle is intended for use within the professional health care facility environment as well as the home health care environment. The shower cradle is not suitable for any other environment or special environments.

The shower cradle is not intended to be used in environments where there are rapid changes in the environmental temperature and humidity during intended use.

2 Components/Key Parts

Please see below to familiarise yourself with the components of the shower cradle. If you have not received all the components listed below contact your local dealer immediately – contact details are provided on the last page of this manual.

Item	Description
1.	Tilt In Space Lever
2.	Removable Mesh Fabric
3.	Stainless Steel Frame
4.	Height Adjustable Braked Castors
5.	Height Adjustable Clips
6.	Gas Struts
7.	User Manual

Table 2




Figure 1

3 Accessories and Spare Parts

The list below shows a list of approved compatible accessories and spare parts that can be ordered at any time during the life of the shower cradle. Contact your local dealer to order these parts. The serial number of your shower cradle can be found on the label on the base of the metalwork.

REF 155003 & 155004 = Small

REF 155005 & 155006 = Medium

REF 155007 & 155008 = Large

3.1 Accessories

The approved accessories are listed below, make sure you know the size of your cradle to ensure the correct parts are ordered. Accessories ordered after purchase will have installation instructions with them.

1.	Adjustable Head Cushion 155805
2.	Adjustable Lap Strap (Small/Medium) 155806
3.	Adjustable Lap Strap (Large) 155807
4.	Commode Tray & Pan Kit (Small/Medium) 155800
5.	Commode Tray & Pan Kit (Large) 155801
6.	Upper Side Rail Kit (Small) (Pair) 155802
7.	Upper Side Rail Kit (Medium/Large) (Pair) 155803
8.	Lower Side Rail Kit (Small/Medium/Large) (Pair) 155804

Table 3

3.2 Spare Parts

The approved spare parts are listed below, make sure you know the size of your cradle to ensure the correct parts are ordered. The correct way to replace the below components is listed in the service manual.

1.	Gas Strut Kit (Pair) 155900
2.	Tilt In Space Hand Lever Kit 155901
3.	Replacement Castor Kit (Single) 155902
4.	Replacement Castor Kit (Set of 4) 155903
5.	Castor Clips (Set of 4) 155904
6.	Small Mesh Cradle Fabric Without Aperture 155140
7.	Small Mesh Cradle Fabric With Aperture 155141
8.	Medium Mesh Cradle Fabric Without Aperture 155142
9.	Medium Mesh Cradle Fabric With Aperture 155143
10.	Large Mesh Cradle Fabric Without Aperture 155144
11.	Large Mesh Cradle Fabric With Aperture 155145

Table 4

4 Assembly Instructions

4.1 Unpacking

The shower cradle will arrive to you in a robust box, please be careful when removing the components from the box. Please read the user manual in full before using this product.



This user manual should be kept safe for future reference.

The shower cradle has been specifically designed to be used in both the professional and home health care environments.

No matter the environment, health and safety factors should be considered to ensure the safety and essential performance of this product and to avoid unnecessary damage or injuries to people within the area.

The environment in which the shower cradle is situated in, whilst carrying out a transfer, is required to be on a flat surface with no steep inclines.



When using a sharp knife, be careful not to damage the product.

This section will summarise the layout of the product packaging and what is included in the box. It is recommended a knife is used for smoother unpacking. The shower cradle is packed into a single palletised box (1,200x800x1,200 mm), weighing approximately 16 kg, excluding the weight of the pallet.

Using a knife to open the box, the internal packaging must be removed to access the Product. It will include all the following components.

1. Fully assembled product
2. User manual
3. Accessories

5 Operating Instructions



This product may be heavy for users and will need to be lifted with care, please operate this product with the use of the castors wherever possible to prevent any injury.

Before use, you should familiarise yourself with the use of the shower cradle, the procedure for tilting the product forwards and backwards and the operation of castor brakes.

Appropriate training in lifting and handling procedures should be undergone by any person operating this product along with any other additional aid equipment, for their own and the client's safety and comfort. These instructions are designed to cover the method of using the product.

5.1 Operating the Tilt Lever



Operation of the tilt lever should only be done when someone is present in the shower cradle. When empty the shower cradle should also be in the fully tilted forward position.



Figure 2



Figure 3



Figure 4

Figure 2 shows the location of the tilt lever. To tilt the shower cradle, place both hands on the rubber handlebar (Figure 3) and squeeze the lever (Figure 4). Whilst the lever is engaged, push down, or pull up on the handle to tilt the shower cradle back or forward, respectively. When the desired tilt position is achieved, release the tilt lever.

5.2 Operating the Braked Castors

The castors on the shower cradle are individually braked using your foot.



Make sure appropriate footwear is worn when using this product. Do not try to push the product when all four castors' brakes are activated.



Figure 5



Figure 6



Figure 7

As per Figure 5, place your foot on the red section of the castor and push down (Figure 6). The red plastic will remain pointing upwards, indicating the brake is now activated. Repeat the process for all castors to ensure that the product is fully secure in position. To release the brakes, press down on the top of the red plastic (Figure 7). The red plastic will remain in the down position. Repeat this for all remaining castors.

5.3 Attaching & Detaching the Shower Cradle Mesh

The mesh can easily be detached and attached to the main frame cradle. This allows for easy cleaning and replacement of the mesh. For detailed instructions on cleaning, washing, and drying of the mesh fabric see the section in this user manual.



The mesh must be fitted correctly to ensure the safe and prolonged use of this product.



Figure 8



Figure 9

Figure 8 shows the metalwork with no mesh fabric attached. To attach the mesh fabric, start at the bottom and slide the pocket over the metalwork (Figure 9), ensuring that the label is visible and facing towards the four castors.



Figure 10



Figure 11




Figure 12

Continue working up the metalwork connecting the buckles one at a time (Figure 10). When you get to the top connect the three 25mm buckles and pull the belts to tighten the mesh (Figure 12).

Once all the buckles are connected, pull the belts to tighten up at each location. The left-over belt should then be placed inside the belt pockets so that they are not hanging down (Figure 11).

If you have purchased an aperture option, then there will be a Velcro section at either side of the aperture section. Pull this tight and connect the Velcro after all buckles have been connected and tightened.

To remove the mesh aperture, loosen all buckles and then unclip them, in order, from the top, down to the bottom.

5.4 Manually Adjusting the Product Height



Never attempt to adjust the height of the product when someone is sitting in the shower cradle. All four castors must be set to the same height before use. Make sure that the plastic retaining clip is fully inserted, the pin must pass through the tube.



Adjusting one castor at a time, remove the plastic retaining clip (Figure 13). Slide the castor bar up or down to the desired height setting ensuring the hole on the tube lines up with the hole on the castor bar (Figure 14). Place the plastic retaining clip back in the hole ensuring that it fully goes through the tube (Figure 15). Repeat this for the three remaining castors.

5.5 Using the Commode Pan



Make sure that the commode pan is empty before attempting to tilt the shower cradle back.



Figure 16



Figure 17



Figure 18

The plastic commode pan self-locates onto the metal commode tray and can easily be removed for emptying the contents and cleaning, by simply sliding the pan out. Figures 16, 17 & 18 show the process of fitting the commode pan to the tray.

5.6 Positioning the Head Cushion

The head cushion is positioned at the top of the cradle with the belt and buckle positioned just above the handle.



Figure 19



Figure 20



Figure 21

Figure 19 shows how to position the cushion. The buckle is then wrapped around the frame and clipped together (Figure 20). Pull the belt to tighten up the cushion and slide the zip tie to prevent the leftover belt from hanging down (Figure 21).

5.7 Accessories



If accessories are bought at the same time as the shower cradle then they will come pre-installed. However, if they are ordered separately then they must be attached by a competent person only. A 5mm Allen key is the only tool required for attaching the accessories and the side rails and commode tray kits are the only accessories that require installation.

For further information contact your local authorised dealer.

6 Technical Specification

6.1 Shower Cradle Dimensions

All dimensions are shown in the table below.



6.2 Specifications

Dimensions – millimetre / (inch)			
Cradle Size	Small	Medium	Large
A – Seat Height	592 mm / (23.3 ")	592 mm / (23.3 ")	592 mm / (23.3 ")
B – Cradle End Height	176 mm / (6.9 ")	116 mm / (4.6 ")	112 mm / (4.4 ")
C – Handle Height	940 mm / (37 ")	1,050 mm / (41.3 ")	1,066 mm / (42 ")
D – Overall Height	1,106 mm / (43.5 ")	1,214 mm / (47.8 ")	1,256 mm / (49.4 ")
E – Base Frame Length	800 mm / (31.5 ")	800 mm / (31.5 ")	800 mm / (31.5 ")
F – Back Length	620 mm / (24.4 ")	750 mm / (29.5 ")	800 mm / (31.5 ")
G – Seat Length	300 mm / (11.8 ")	350 mm / (13.8 ")	400 mm / (15.7 ")
H – Leg Length	360 mm / (14.2 ")	450 mm / (17.7 ")	460 mm / (18.1 ")
I – End Length	150 mm / (5.9 ")	150 mm / (5.9 ")	150 mm / (5.9 ")
J – Internal Width	440 mm / (17.3 ")	440 mm / (17.3 ")	520 mm / (20.5 ")
K – Overall Width	502 mm / (19.8 ")	502 mm / (19.8 ")	582 mm / (22.9 ")
Upright Length	950 mm / (37.5 ")	1,134 mm / (44.6 ")	1,134 mm / (44.6 ")
Reclined Length	1,260 mm / (49.6 ")	1,508 mm / (59.4 ")	1,508 mm / (59.4 ")
Turning Circle	836 mm / (32.9 ")	836 mm / (32.9 ")	886 mm / (34.9 ")
Castor Size	125 mm / (4.9 ")	125 mm / (4.9 ")	125 mm / (4.9 ")

Table 5

Weights – kilogram / (pounds) / (stone)	
Safe Working Load (SWL)	160 kg / (352 lb) / (25 st)
Product & Packaging	16 kg / (35.2 lb)
Product	15 kg / (33 lb)

Table 6

6.3 Expected Product Lifetime

7 years, depending on usage and compliance to maintenance, servicing, and inspections.

7 Environmental - Storage and Operating Conditions

The shower cradle is intended for internal use within normal environmental conditions.

It is not intended to be used in environments where there are rapid changes in the environmental temperature during intended use.

7.1 Normal Operating Conditions

+5°C to +40°C (41°F to 104°F) at a relative humidity between 15% to 90% RH, non-condensing but not requiring a water vapour pressure greater than 50hPa and atmospheric pressure between 700hPa to 1060hPa

7.2 Shipping and Storage Conditions

-25°C to +5°C (-13°F to 41°F) with any humidity level.

+5°C to +35°C (41°F to 95°F) at a relative humidity up to 90%.

+35°C to 70°C non-condensing at a water vapour pressure up to 50hPa

8 Disposal

When the shower cradle has completed its life cycle and can no longer perform to its intended use safely, the product must be decommissioned by an approved Service Engineer. The following specifies the importance of correct disposal procedure including local laws and being environmentally friendly.

Please observe the local laws on recycling and respect the current laws for disposal within the community the product is being used within. If there is any uncertainty of the below guidelines, contact your local authorities to determine the proper method of disposal of potentially biohazardous parts and accessories.

The relevant components utilised in the manufacture of the product that can be recycled at the end of the product life are:

Fully Recyclables	Consideration when Recycling
Steel frame and accessories	Actuators
Castors	Rubber handlebar covers
Initial packaging of the product (cardboard)	Plastic fittings
Metallic fixings – screws etc.	
Mesh fabric	

Table 7



The product may be contaminated and must be disinfected before recycling or disposal. See section 'Cleaning' for details of how to do this.

9 Fault Finding

Should a problem arise with the use of the shower cradle, review the Table below.

Find the fault and complete the recommended solution.

If the fault is not found and/or the solution does not correct the problem, contact your local authorised dealer immediately – contact details are provided on the last page of this manual.

Fault	Action
Castors do not operate	Ensure the brake has been fully disengaged, Check that the castors are free running and that they are touching the ground. If this does not correct the fault, then contact your local authorised dealer immediately so the product can be checked to ensure proper continued operation.
The shower cradle doesn't tilt	Ensure the lever touches the handlebar when pressed fully. Check that the cables are in good condition and there are no kinks in the cable. If this does not correct the fault, then contact your local authorised dealer immediately so the product can be checked to ensure proper continued operation.
The mesh fabric is loose	Check that every buckle has been connected then pull the belt on each buckle to tighten at every point. If this does not correct the fault, then contact your local authorised dealer immediately so the product can be checked to ensure proper continued operation
Movement is compromised by excessive force required to move	Check the castors are not fouled with debris that would restrict its movement (clean out, test and return to service) If this does not correct the fault, then contact your local authorised dealer immediately so the product can be checked to ensure proper continued operation.
Brakes do not operate	Ensure that the brake mechanism is fully engaged, check to see that all castors are touching the ground. If this does not correct the fault, then contact your local authorised dealer immediately so the product can be checked to ensure proper continued operation.

Table 8

10 General Inspection and Maintenance

10.1 Service



No service is to be carried out on the shower cradle while transferring a person to reduce the risk of injury.

Service must be completed by a Prism Medical UK authorised Service Engineer only.

Do not attempt to service the product yourself, or warranty is void.

To ensure the safety and continued good function of your shower cradle, routine service must be performed. Service should be completed by an approved service engineer every one year to ensure the products required standard is maintained. The service history of the product should be documented each service in the Service Log at the back of this User Manual.



When the product is serviced, the service checklist must be completed.

Service Manual Document Number: 995147.

Contact your local authorised dealer if you:

- Need more information.
- Have any questions about the use or service of your shower cradle.
- Notice any change in the performance.
- Want to report an unexpected occurrence.
- Want to arrange a service.
- Need to ascertain necessary information for replacement parts and components.

Contact details of your local Prism Medical dealer are shown on the last page of this manual

10.2 Inspection

Inspection is to be completed prior to each use by the user of the shower cradle.



Should any of the components in the table below fail the inspection, DO NOT use the shower cradle. Contact your local authorised dealer for service – contact details are on the last page of this manual.

Ensure all component inspections in the Table below are completed prior to each use of the shower cradle.

Check List before Use

Component	Service/Inspection required
Generic	Visual inspection of the external of the shower cradle. Significant damage that may affect the function of the shower cradle along with a clear safety hazard is unacceptable.
	Check the Labelling on the metalwork and mesh to ensure they are all still legible, this includes the Serial Number and other important markings. If labels are not legible, then contact your local authorised dealer immediately.
	Check all main nuts and bolts to see if they are loose, if they are not tight or you have concerns, then contact your local authorised dealer immediately.
	Examine the tilt lever and test that it still operates smoothly.
Mesh Fabric	Inspect the mesh fabric for any damage including cuts and breaks. DO NOT use if there is fraying or damage to any part of the fabric.
Brakes	Ensure the brakes function when activated.
	Ensure the product can move freely when the brakes are deactivated. (Ensure they aren't rubbing on the wheels).
Castors and Wheels	Examine the products' castors for signs of damage or wear, ensure they run freely and are not cracked.
	Check to see all castors are flat on the floor and shower cradle is stable.
	Ensure the product can rotate freely with ease.

Table 9

11 Cleaning

Please follow the cleaning guidelines below on cleaning and disinfecting the shower cradle.

11.1 Washing & Drying of the Mesh Fabric



Follow the instructions from section 5 to remove the mesh fabric from the cradle.

The mesh fabric can then be machine washed using a non-biological soap powder up to a temperature of 80°. Prism Medical recommends air drying the fabric only.

11.2 General Cleaning



It is recommended to clean the shower cradle and accessories before use by a different person, reducing the risk of cross-contamination.

The exterior of the shower cradle can be cleaned using a damp soapy cloth for general cleaning duties. Please ensure the cloth is damp and not wet. Ensure the exterior of the device is dry after cleaning. Dry using a clean dry cloth.

11.3 Disinfecting (if necessary)

Should the shower cradle require a more thorough clean, the use of the Actichlor™ disinfectant product (which is widely available in tablet form and used throughout the health care industry) is recommended.



Follow the manufacturer's safety instructions for the use of the cleaning product before use to ensure safe use for the operator and the patient.
Ensure the cloth is damp before the cleaning process.

Application is through a clean damp cloth applied to wipe the product down. Use in the following dilutions to ensure an effective clean:

- Actichlor™ dissolvable chlorine tablets provide a concentration of 1000 ppm of available chlorine (0.1%) per 1 tablet
- 1 tablet (1.7g formed tablet (x1) will create a virucidal solution, diluted in 1 litre of water to provide effective means to clean a "dirty" product. This is also ideal for use after an outbreak of the Norovirus/winter vomiting and can be used as a precaution against C.Diff. It is effective against viruses, bacteria, spores, yeasts and moulds.
- The contact time against the outer components of the product should be for 5 minutes to prevent any virucidal infections without a degradation to the functionality of the product. 5 minutes is a recommended contact time. The product can withstand a longer contact period but the 5-minute recommendation as a minimum must be followed to provide an effective cleaning regime.
- Blood spills should be dealt with by an increased concentration of the solution – please refer to the instructions on the manufacturer's product labelling.

Dilution chart					
Product used as	Product condition	Concentration (ppm)	Dilution qty* (l)	Tablets per 1l (0.26gal)	Contact time (minutes)
Bactericidal	Clean	200	5 (1.32gal)	1	1
	Dirty	1000	1 (0.26gal)	1	5
Yeasticidal	Clean	200	5 (1.32gal)	1	1
	Dirty	1000	1 (0.26gal)	1	5
Fungicidal	Clean	2000	1 (0.26gal)	2	15
	Dirty	5000	1 (0.26gal)	5	15
Mycrobactericidal	Clean	1000	1 (0.26gal)	1	15
	Dirty	5000	1 (0.26gal)	5	15
Virucidal	Clean	500	2 (0.53gal)	1	5
	Dirty	1000	1 (0.26gal)	1	5
Sporcidal (C.Diff)	Clean	1000	1 (0.26gal)	1	10
	-	-	-	-	-
Sporcidal	Clean	5000	1 (0.26gal)	5	10
	-	-	-	-	-

* Dilution is made with water. DO NOT dilute within any other medium.

- When diluted in water, one tablet gives 1000ppm of available chlorine.
- The concentration of the solution depends upon whether the object being cleaned is noticeably dirty (indicated in the table by "Product condition").

Table 10

Handling and storage safety precautions when using this cleaning agent:

Advice on Safe Handling



- Avoid contact with skin and eyes.
- Do not breathe dust/fumes/gas/mist/vapours/spray.
- Use only with adequate ventilation.
- Wash hands thoroughly after handling.
- Mixing this product with acid or ammonia releases chlorine gas.

Hygiene Measures

Handle in accordance with good industrial hygiene and safety practice. Remove and wash contaminated clothing before re-use. Wash face, hands, and any exposed skin thoroughly after handling.

Conditions for Safe Storage



- Keep out of reach of children.
- Keep container tightly closed.
- Store in suitable labelled containers.
- Storage temperature: 0-25°C (32-77°F).

Individual Protective Measures

Hand protection: Gloves

Dissolve

Dissolve in cold water – With no agitation, 1 tablet will take approximately 10 minutes to fully dissolve in the water used.

The information above has been extracted from the Actichlor™ MSDS (Manufacturers Safety Data Sheet). For a full review of the data please follow the link below:

<http://www.nhsggc.org.uk/media/236215/msds-actichlor-plus.pdf>

12 Warranty

This guarantee does not affect or in any way limit your Statutory Rights.

1. The manufacturer guarantees the Freeway Flow Shower Cradle, against failure within the period of 60 months from the date of purchase by virtue of defects in material or workmanship. This is excluding the shower cradle mesh which has a warranty of 12 months from the date of purchase.
2. The liability of the manufacturer under terms of this guarantee shall be limited to the replacement or the defective part(s) to the sales distributor, dealer, agent, person, or entity which purchased the equipment from the manufacturer. In no event shall the manufacturer incur liability for any consequential or unforeseeable losses.
3. This equipment guarantee shall be void if the equipment is not serviced by Prism Medical UK or its authorised service engineers, in accordance with manufacturer's recommendations, or if any unauthorised persons carry out work on the equipment.
4. This guarantee does not apply to failure attributable to normal wear and tear, damage by natural forces, user neglect or misuse or to deliberate destruction.
5. Do not attempt to service the product yourself, or warranty is void.

13 Service Record History

Complete this section after each service, repair inspection and/or maintenance.

Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	

Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	

Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	
Date:	Time:
Service Type: Service Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Other <input type="checkbox"/>	
Completed By: (printed name) (signature)	
Company:	
Remarks & Actions Taken:	
Product Left in A Safe & Usable Condition: Yes <input type="checkbox"/> No <input type="checkbox"/> (if no explain in actions above)	

User notes:

Dealer / service contact details:

Manufacturer contact details:

Prism Medical UK

Unit 1, Tir Llwyd Industrial Estate,
St Asaph Avenue, Kinmel Bay,
Conwy, United Kingdom, LL18 5JZ

Telephone Number: 01924 840100

Disclaimer

While every effort has been made to ensure the accuracy of information contained in this user manual, no liability can be accepted by Prism Medical UK for any errors or omissions. Prism Medical UK operates a policy of continuous improvement. Specifications and other data are subject to change without notice.

www.prismmedical.co.uk

