

Falls Recovery Sling



CONTENTS

1.0	Safety Instructions and Warnings	3
1.1	Introduction	3
1.2	Manufacture	3
1.3	European Authorised Representative.....	3
1.4	Symbols Used.....	4
1.5	Contraindications / Limitations.....	5
1.6	Intended Use.....	5
1.7	Warning Notes	6
1.8	Operating Environment	6
2.0	Components/Key Parts	7
3.0	Assembly Instructions	8
3.1	Unpacking	8
3.2	Unboxing.....	8
4.0	User Instructions.....	9
4.1	Initial Safety Checks	9
4.2	Positioning the Person onto the Sling.....	9
4.3	Preparing the Sling to Move the Patient along the Floor	10
4.4	Preparing the Sling for Hoisting	10
4.5	Hoisting the Person into the Chair.....	12
5.0	Technical Specification.....	13
5.1	Sling Dimensions	13
5.2	Specifications	14
5.3	Expected Product Lifetime	14
5.4	Standards Applied	14
6.0	Environmental - Storage and Operating Conditions.....	15
6.1	Normal Operating Conditions	15
6.2	Shipping and Storage Conditions	15
7.0	Disposal.....	15
8.0	General Inspection and Maintenance.....	16
8.1	Service.....	16
8.2	Inspection.....	17
9.0	Cleaning.....	18
9.1	General Cleaning.....	18
9.2	Disinfecting (if necessary)	18
9.3	Drying.....	18
10.0	Warranty	19
11.0	Inspection Record History.....	20

1.0 Safety Instructions and Warnings

1.1 Introduction



Please read and understand this manual in its entirety before using your Freeway Falls Recovery Sling. This manual includes all variants of the Freeway Recovery Falls Sling, this includes:
Freeway Recovery Falls Sling – Loop Attachment
Freeway Recover Falls Sling – Clip Attachment

The information in this manual is important for the safety of anyone near the Falls Recovery Sling and must be read and understood to help prevent injuries. It is also crucial to the proper operation and maintenance of the Sling.

Should any questions arise from reviewing this manual, contact your local authorised representative.

Failure to comply with warnings in this manual may result in; injury to the operator and/or client and/or damage to the Sling or related components.

If, during the use of this sling or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

Store this manual with the documents included with the Sling. Contents of this manual are subject to change without prior written notice.



Do not attempt to use this equipment without first understanding the contents of this manual.



Unauthorised modifications on any Prism Medical UK product may affect its safety. The manufacturer will not be held responsible for any accident, incident or deficiencies of performance that occur as a result of any unauthorised modification to its products.

1.2 Manufacture

The Slings is manufactured at the address below:



Prism Medical UK
Unit 1, Tir Llwyd Industrial Estate, St Asaph Avenue, Kinmel Bay, Conwy, LL18 5JZ
Telephone number: 01924 840 100

1.3 European Authorised Representative

The address of the European Authorised Representative for this product:



European Healthcare & Device Solutions (Ireland) Ltd.
Stratton House, Bishopstown Road,
Cork, Ireland.
T12 Y9TC.
Telephone number: +353(86)2280846

1.4 Symbols Used

The Table below includes all Symbols from BS EN ISO 15223-1:2016 that can be found in this manual and on the product and what they represent. Refer back to this table when you are unsure of what a symbol represents.





















	Consult instructions before use		Caution – see instructions for use
	Manufacturer	SWL	Safe Working Load represents the maximum load rated for safe operation
	Packaging indicator – This way up		Date of manufacture
	Serial number		Packaging indicator – Keep dry
	Please observe local laws on recycling		For internal use only
	Temperature range	IP_{N₁}N₂	Degree of protection provided by enclosure. N ₁ : Ingress of particles N ₂ : Ingress of water
	Humidity range		Medical Device
	Catalogue number		European Authorised Representative
	Machine wash hot		Do not bleach
	Do not Iron		Tumble Dry Low Heat
	Loop style attachment for sling		Clip Style attachment for sling

Table 1

1.5 Contraindications / Limitations

There are no known “contraindications” associated with the usage of the Freeway Falls Recovery Sling, provided they are used as per manufacturer’s recommendations and guidelines. However, it is recommended that a client specific assessment is completed by a trained and knowledgeable health care professional to determine the method of use.

Prism Medical UK does not recommend a required number of care givers for the use of our products. This information and recommendation can only be provided after a thorough personalised, case specific assessment, as there are many factors that can influence these decisions.

1.6 Intended Use



For internal use only.

The Freeway falls recovery sling has been developed to allow carers to safely transfer a person that needs assistance immediately after a fall. The sling will be generally used with a mobile person who has happened to have had a fall. It is important that you have assessed the person to see that they haven’t sustained any injuries and that medical attention isn’t required. The sling is designed to fit all sizes and shapes of user, up to the safe working load of the sling. The sling is designed for ages 4 and above.

The Freeway Falls Recovery Sling allows for the fallen person to be rolled onto the sling in a lying position. Once onto the sling the fallen person can be slid into a more open area if required, for the suitable positioning of the hoist and user.



Please note: The Yellow side of the sling is the sliding side, which is the underside. This is a low friction material and is fitted to the sling to aid the sliding of the person in the sling. The yellow straps and yellow attachment loops on the sling are used for the sliding of the person on the sling only. These are not used for the hoisting of the person.

Once in a suitable position within the room, a mobile hoist, ceiling track hoist or any other suitable hoisting equipment can be introduced and the fallen person can be hoisted up off the floor safely back into a chair or directly back onto a bed. The carer(s) must be assessed to be a competent person, who must be trained in the suitable application and fitting of the slings. A video showing the safe use of the sling can be found at www.prismmedical.co.uk.



The Freeway Falls Recovery Sling has a SWL (safe working load) of 200 kg, this load must not be exceeded. When using this sling with other medical equipment, the sling SWL must not be exceeded when used with a hoist with a higher specified load. sling is used under instruction, and the operation of the sling is undertaken by a trained carer.

A risk assessment must be performed before using any other manufactured sling, to ensure safe use can be established. A risk assessment must also be carried out when using this sling with other medical devices.

The carer is required to undertake basic movement and handling training to ensure they are competent in using the sling for its intended use.

Summary:

- Assists in moving a fallen person by sliding them into a more adequate lifting position and allows them to be lifted by fixing onto a medical hoist.
- Helps to transfer a fallen person from one destination to another (i.e. bedroom to bathroom)
- Maximum user weight: 200 kg
- Low friction material to allow the person to be slid across the floor.
- Padded sling ensures comfort at all times for the user.
- Multiple loops/clips available for hoist attachment, choose the best suited attachment for the situation.
- Yellow straps are used for sliding the person only, not to be used for lifting.

1.7 Warning Notes



Your sling has been manufactured and tested to exceed BS EN ISO 10535:2006 – Hoists for the transfer of disabled persons. This does not mean that it can be used without care. All operators should have read the operating instructions and appreciate this warning section.

- The sling is intended to be used on a flat or levelled surface when used in conjunction with a mobile hoist.
- Under no circumstance should the sling be put in control of a person who has not been properly trained in the use and care of this equipment. Failure to adhere to this warning may result in serious injury to the operator and/or the individual.
- In facilities where more than one operator will be responsible for using the sling, it is imperative that all such members be trained on the sling prior to use. A training program should be established by the facility to acquaint new operators with this equipment.
- Your guarantee is void if persons unauthorised by Prism Medical UK perform work or modification to the sling.
- To maintain optimum function and safety, the sling should be inspected and maintained on a regular basis. See section 'General Inspection, Maintenance and Cleaning' within this user manual.
- In areas where children are prone to be present, be vigilant when using the sling.
- Any accessories used with the sling, such as a mobile hoist, should be checked to ensure that they are in good working order. Check for signs of wear to each component prior to use. Report any unusual wear to your local authorized dealer.
- The sling is certified to a maximum load of 200kg. Do not exceed the maximum rated load.
- Ensure that a clear space is maintained around the sling before performing a transfer, check for and move all obstacles out of the way.
- Protecting the people present, when using the sling with a medical hoist, visually monitor the clipped/looped connection points during transfer stages, so the sling remains firmly attached to the connection pins/hooks.
- To reduce the risk of unintended use, when the sling is not in use remove the sling from any medical hoist to prevent entrapment or strangulation should the device be tampered with. The sling can be placed back in the
- Between the sling and any medical hoist, the lowest maximum load shall always be used.



You may need to seek specialist advice on how to assist some people with specific moving and handling needs. Sources of advice include, but is not limited to, professional bodies and organisations, occupational therapist, physiotherapists, manual handling advisers and ergonomist with experience in health and social care.

1.8 Operating Environment

The sling is intended for use within the professional health care facility environment as well as the home health care environment. The sling is not suitable for any other environment or special environments.

The sling is not intended to be used in environments where there are rapid changes in the environmental temperature and humidity during intended use.

2.0 Components/Key Parts

Please see below to familiarise yourself with the components of the Falls Recovery Sling. The images below show the contents of the sling. If you have not received all the components contact your Prism Medical UK dealer immediately – contact details are provided on the last page of this manual. The images below show the two types of sling available, this includes the clipped type attachment and the loop type attachment

Item	Description
1	Main sling body – patient sits here
2	Low friction material – for sliding
3	Sling main label
4	Yellow straps – for sliding
5	Lifting Straps – for lifting
6	Lifting Clips – for lifting
7	Sling head support
8	Non-slip bag
9	Yellow Pull Straps – for sliding

Table 2

freeway

PLEASE REFER TO USER INSTRUCTIONS
EXAMINE SLING BEFORE USE. DO NOT USE IF DAMAGED

THIS SLING MAY BE USED WITH SUITABLE HOISTS TO A MAX SWL OF:
200kg/440lb

Prism Medical UK Ltd
Unit 1, Tir Lwyd Industrial Estate
St Asaph Avenue, Kinnel Bay, Rhyl
Conwy LL18 5JZ

European Healthcare & Device Solutions (Ireland) Ltd
Stratton House, Bishopstown Rd
Cork T12 Y9TC

+44 (0)1924 840100 | www.prismmedical.co.uk

LOLER EXAMINATION RECORD

YEAR	ISSUE DATE:											
1	1	2	3	4	5	6	7	8	9	10	11	12
2	1	2	3	4	5	6	7	8	9	10	11	12
3	1	2	3	4	5	6	7	8	9	10	11	12
4	1	2	3	4	5	6	7	8	9	10	11	12
5	1	2	3	4	5	6	7	8	9	10	11	12
	1	2	3	4	5	6	7	8	9	10	11	12

Figure 1

Image of Looped Style Sling




Figure 2

Image of Clipped Style Sling



Figure 3



Figure 4

3.0 Assembly Instructions

3.1 Unpacking



The Falls Recovery Sling could arrive to you in a robust box, inside the box the Falls Recovery Sing has its own non-slip bag, do not discard the bag, this bag is used to aid when seating a person onto a chair during recovery, as well as being a bag to store the sling when not in use.

This user manual should be kept safe for future reference.

The sling has been specifically designed to be used in both the professional and home health care environments.

No matter the environment, health and safety factors should be considered to ensure the safety and essential performance of the sling and to avoid unnecessary damage or injuries to people within the area of the sling.

3.2 Unboxing



When using a sharp knife, be careful not to damage the product.

This section will summarise the layout of the sling packaging and what is included in the Box. It is recommended a knife is used for smoother unpacking of the sling. Usually, (if two slings have been purchased), they will be packed into a single box (400x475x610), weighing approximately 8kg.

4.0 User Instructions



Before use, you should familiarise yourself with use of the sling, the procedure for placing a person in the sling, securing them in place correctly, and attaching the sling onto a medical hoist.

Appropriate training in lifting and handling procedures should be undergone by any person using a sling, for their own and the client's safety and comfort. These instructions are designed to cover the method of using the sling.

See below for details on how to use the sling for its intended purpose:

The recovery falls sling should be used when a person has been found lying on the floor.

Images below refer to the looped version of the Falls Recovery Sling, but all images can be used as reference for both sling types unless stated otherwise.

4.1 Initial Safety Checks

First of all you should:

- Check the person for any serious injuries.
- If the person is seriously injured, the emergency services must be contacted.
- If the person is slightly injured, perform basic first aid and make the person comfortable.
- Assess the surroundings of the person, ensure that the surrounding environment is safe.
- Ensure that the person is in a suitable condition that they can be moved and/or hoisted.
- Summon assistance to move the person onto the falls recovery sling.

4.2 Positioning the Person onto the Sling

Once the steps in 4.1 have been complete, you can now move the person onto the falls recovery sling.



It is required that two carers move the patient onto the falls recovery sling.

- Using the correct manual handling techniques, carefully roll the person onto their side as shown. (Figure 5)
- Place the sling onto the floor with the yellow low friction material face down.
- Double fold half of the sling back on top of itself (only the black material should be on show) and slide the sling up towards the back of the patient. (Figure 6)
- Roll the patient onto their other side, with the sling now directly underneath them.
- You should now be able to pull out the folded part of the sling from underneath the patient.
- The patient should be fairly central within the sling. (Figure 7)



Figure 5




Figure 6




Figure 7

4.3 Preparing the Sling to Move the Patient along the Floor

Once the steps in 4.2 have been complete, you can now set up the sling for patient transfer along the floor to an adequate position for hoisting.

- Attach the pull straps (found separately in the bag, these are used to slide the person in the sling) onto the small yellow hook straps, this can be done in either of the two ways shown, whichever is preferred. (Figure 8 and Figure 9)
- Fasten the red buckle as shown in figure 10. This will provide extra security when moving the patient along the floor.
- Once the straps are all set, the carer can now slide the patient along the floor to a suitable location for hoisting, or location to await emergency services if required. (Figure 11)
- Always be vigilant of any obstacles while moving the patient. Move them slowly and carefully to ensure they remain comfortable and safe.



Figure 8



Figure 9




Figure 10



Figure 11

4.4 Preparing the Sling for Hoisting

Once the steps in 4.3 have been complete, you can now set up the sling for hoisting the patient.

- When the sling is ready to be hoisted, the yellow pull straps can be removed, these should be kept in a safe place until the sling is no longer in use. Once finished with the sling, remember to place the yellow straps back into the bag for safe keeping.
- Locate the centre straps on the sling and feed them between the thighs of the patient.
- Once the two centre straps are in hand, the loops/clips can be fed through the hoisting straps on their respective sides. See figure 12 for reference.



Figure 12

- The four main hoisting straps can now be fitted onto the hoist attachment point.
- For the looped version of the Recovery Falls Sling, loop the straps onto the hoist fixing point as shown below:



Figure 13



Figure 14



Figure 15

- For the clipped version of the Recovery Falls Sling, clip the clips onto the hoist fixing point as shown below:



Figure 16



Figure 17



Figure 18

- For further information on sling attachment to the hoist, refer to the user manual of the hoist used, that document should show how to attach a sling to that specific product.

4.5 Hoisting the Person into the Chair

Once the steps in 4.4 have been complete, you can now move the patient in the hoist and set them down into a chair. A similar process can be followed when hoisting the person into a bed or sofa.

- Raise the lifting arm of the hoist to raise the patient of the ground within the sling. Always double check that the patient is secure within the sling. (Figure 19 and Figure 20)
- The non-slip bag, which the sling is provided within can now be used to assist in placing the person into the chair. (Figure 21)
- Place the strap of the bag over the back of the chair and place the bag onto the seat. (Figure 21)
- Transfer the person within the hoist directly above the chair.
- Rotate the person within the sling to position them into the chair during lowering. (Figure 22)
- Lower the person into the chair and disconnect the sling from the hoist. The person is now safely positioned back into their chair. (Figure 23)



Figure 19



Figure 20




Figure 21



Figure 22



Figure 23

5.0 Technical Specification

5.1 Sling Dimensions

All dimensions below are shown in millimetres.

Looped Version




Figure 24

Clipped Version



Figure 25

5.2 Specifications

Dimensions	Clipped Version Sling	Looped Version Sling
Maximum Length (Including Head Support)	1840mm	1480mm
Maximum Width	1180mm	1100mm
Thickness	35mm	35mm

Table 3

Weights		
Safe Working Load (SWL)	200 kg	200 kg
Full Product Weight	1.7 kg	1.7 kg
Bag Weight (with Straps)	1.3 kg	1.3 kg

Table 4

Other Specs		
Material Type	Super Soft Spacer	Super Soft Spacer
Sizes Available	One Size Fits All	One Size Fits All
Head Support	Yes	Yes
Bag Available	Yes	Yes

Table 5

5.3 Expected Product Lifetime

5 years, when manufactures use instructions and laundering requirements are adhered to.

5.4 Standards Applied

The standards that have been applied to the device are as follows:

- BS EN ISO 10535:2006 – Hoists for the transfer of disabled persons – Requirements and test methods.
- EN ISO 15223-1:2016 Medical device Symbols
- BS ISO 20417:2021 Information supplied by the manufacture

6.0 Environmental - Storage and Operating Conditions



The sling is intended for internal use within normal environmental conditions.

It is not intended to be used in environments where there are rapid changes in the environmental temperature and humidity during intended use.

6.1 Normal Operating Conditions

+5°C to +40°C (41°F to 104°F) at a relative humidity between 15% to 90% RH, non-condensing but not requiring a water vapour pressure greater than 50hPa and atmospheric pressure between 700hPa to 1060hPa

6.2 Shipping and Storage Conditions

-25°C to +5°C (-13°F to 41°F) with any humidity level.

+5°C to +35°C (41°F to 95°F) at a relative humidity up to 90%.

+35°C to 70°C non-condensing at a water vapour pressure up to 50hPa.

7.0 Disposal



When the sling has completed its life cycle and can no longer perform to its intended use safely the sling must be decommissioned by an approved Service Engineer. The following specifies the importance of correct disposal procedure including local laws and being environmentally friendly.

Please observe the local laws on recycling and respect the current laws for disposal within the community the device is being used within. If there is any uncertainty of the below guidelines, contact your local authorities to determine the proper method of disposal of potentially biohazardous parts and accessories.

The relevant components utilised in the manufacture of the device that can be recycled at the end of the device life are:

Fully recyclables:	Consideration when Recycling:
Plastic Clips	Sling Body
Initial packaging of the device (cardboard)	Sling Bag

Table 6



The product may be contaminated and has to be disinfected before recycling or disposal. See section on 'Cleaning' for details of how to do this.

8.0 General Inspection and Maintenance

8.1 Service



No service is to be carried out on the sling while transferring a person to reduce the risk of injury. Service must be completed by an authorised Service Engineer. Do not attempt to service the product yourself, or the warranty will be void.

To ensure the safety and continued good function of your sling, routine service must be performed on your service. Service should be completed by an approved service engineer every 6 months to ensure the products required standard is maintained. The service history of the product should be documented each service in the service log at the back of this user manual.



A service/inspection must be completed every 6 months after installation of the sling as per the latest LOLER Regulations. The sling has an expected Service Life of 5 Years.

Contact your local authorised dealer if you:

- Need more information.
- Have any questions about the use or service of your sling.
- Want to report an unexpected occurrence.
- Want to arrange an LOLER inspection.

Contact details of your local authorised dealer are shown on the last page of this manual. If your local dealers details aren't available on the last page. Contact the stated manufacture which is also found on the last page of this manual.

8.2 Inspection

Inspection is to be completed prior to each use and after each wash by the user of the sling.



Should any of the components in the table below fail the inspection, DO NOT use the sling. Contact your local authorized dealer for service – contact details are on the last page of this manual.

Ensure all component inspections in the table below are completed prior to each use of the sling.

Check List before Use

Component	Service/Inspection required
Sling Body	Conduct a visual inspection of the sling. Significant damage that may affect the function of the sling along with a clear safety hazard is unacceptable. Damage includes: <ul style="list-style-type: none"> - Cuts - Frays - Tears - Burns - Stitching coming undone - Stretched fabric - Holes in the fabric - Discolouring or bleaching on the material - Significant shrinkage
Straps	Inspect the straps for any significant damage that may affect the function of the sling along with any clear safety hazards. Damage includes: <ul style="list-style-type: none"> - Cuts - Frays - Tears - Burns - Stitching coming undone - Stretched fabric - Holes in the fabric - Discolouring or bleaching on the material - Significant shrinkage
Labels	Check the Labelling on the sling to ensure they are all still legible, this includes the Serial Number and other important markings. If labels are not legible, then it is recommended that the product is decommissioned.

Table 7

9.0 CLEANING

Sensitive laundering of slings will prolong the extent of the expected life cycle of the product. Conversely, improper laundering will reduce the expected life cycle and possibly damage the integrity of the sling, resulting in it having to be discarded.

Please follow the guidelines below on cleaning and disinfecting the sling and the sling bag.

9.1 General Cleaning



It is recommended to clean the Freeway Falls Recovery Sing before it is used by a different person, reducing the risk of cross-contamination.

The sling should be laundered with soap or detergent at 50°C on a standard domestic cycle.

It is recommended that non-biological detergents are used to clean a sling.

Any clips, hooks and buckles should be closed to protect them from damage during a wash cycle.

Avoid using bleach when washing, particularly hypochlorite as it can cause structural damage to the sling. It is also important to use cleaning products with washing powder only, do not include any fabric conditioning as this can cause damage to the sling. Do not overheat the sling during a cycle wash. Temperatures above 80°C will cause shrinkage and structural damage, it is always best to maintain a 50°C cycle wash where possible, but the temperature will need to be increased during disinfection.

It is recommended that the sling is washed after every use to avoid cross-contamination between users.

9.2 Disinfecting (if necessary)

Should the sling require a more thorough clean, the temperature guidelines below should be followed as below:

- 10 minute wash at 65°C or
- 3 minute wash at 72°C

9.3 Drying

We recommend that the slings should air dried within approximately 4 hours within a warm room, this will prolong the service life. Alternatively we recommend the sling can be tumble tried at 50°C if required.



Do not place a sling over a radiator or convectors.

If the cleaning procedures above are followed correctly, along with the correct maintenance and servicing of the product. The sling should achieve its full life expectancy.

10.0 Warranty

This guarantee does not affect or in any way limit your Statutory Rights.

1. Prism Healthcare guarantees the Freeway Falls Recovery Sling, supplied as new, against failure within the period of 12 months from the date of purchase by virtue of defects in material or workmanship.
2. The liability of Prism Healthcare under terms of this guarantee shall be limited to the replacement or the product to the sales distributor, dealer, agent, person or entity which purchased the equipment from Prism Healthcare. In no event shall Prism Healthcare incur liability for any consequential or unforeseeable losses.
3. This equipment guarantee shall be void if the product is not maintained by Prism Healthcare its authorized agents, in accordance with manufacturer's recommendations, or if any unauthorized persons carry out work on the equipment.
4. This guarantee does not apply to failure attributable to normal wear and tear, damage by natural forces, user neglect or misuse or to deliberate destruction.
5. Do not attempt to modify or repair the product yourself, or warranty is void.

11.0 Inspection Record History

Complete this section after each inspection of the product.

Serial Number	
----------------------	--

Date: _____ Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other
Completed by: (printed name). (signature)
Company:
Remarks & Action Taken:
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)

Date: _____ Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other
Completed by: (printed name). (signature)
Company:
Remarks & Action Taken:
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)

Date: _____ Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other
Completed by: (printed name). (signature)
Company:
Remarks & Action Taken:
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)

Date: _____	Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other	
Completed by: (printed name). (signature)	
Company:	
Remarks & Action Taken:	
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)	

Date: _____	Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other	
Completed by: (printed name). (signature)	
Company:	
Remarks & Action Taken:	
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)	

Date: _____	Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other	
Completed by: (printed name). (signature)	
Company:	
Remarks & Action Taken:	
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)	

Date: _____	Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other	
Completed by: (printed name). (signature)	
Company:	
Remarks & Action Taken:	
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)	

Date: _____	Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other	
Completed by: (printed name). (signature)	
Company:	
Remarks & Action Taken:	
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)	

Date: _____	Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other	
Completed by: (printed name). (signature)	
Company:	
Remarks & Action Taken:	
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)	

Date: _____	Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other	
Completed by: (printed name). (signature)	
Company:	
Remarks & Action Taken:	
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)	

Date: _____	Time: _____
Service Type: <input type="checkbox"/> Periodic inspection <input type="checkbox"/> Monthly inspection <input type="checkbox"/> 6-month inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly inspection <input type="checkbox"/> Other	
Completed by: (printed name). (signature)	
Company:	
Remarks & Action Taken:	
Device left in a safe usable condition: YES. <input type="checkbox"/> NO <input type="checkbox"/> (if "NO" explain in remarks the action taken)	

User notes:

Dealer/service contact details:

Manufacturer contact details:

Prism Medical UK

Address: Unit 1 • Tir Llwyd Industrial Estate • St Asaph Avenue • Kinmel Bay • Conwy • LL18 5JZ

Telephone Number: 01924 840 100

Disclaimer

While every effort has been made to ensure the accuracy of information contained in this user manual, no liability can be accepted by Prism Healthcare for any errors or omissions. Prism Healthcare operates a policy of continuous improvement. Specifications and other data are subject to change without notice.

