


mackworth

SA180 STAND AID



SPARE PARTS MANUAL

1. Contents



Description	Page(s)
Spare Parts Summary	3
Part Code Detail	4-9
Notes	10-11
Contacts	12





2. Spare Parts Summary

Part Code	Description
107908	Leg Rod Assembly
119105	Leg Actuator Lead (750mm)
123105	4 Way Handset (Single Battery LED)
123108	4 Way Control Box (Standard)
107904	LH Leg Assembly
107903	RH Leg Assembly
000213	Mast Fixing Screw (Small)
000746	Mast Fixing Screw (Large)
107212	SA180 Base Assembly
107915	Lift Actuator Kit
107928	Leg Actuator Kit
107921	SA180 Boom Assembly
107922	SA180 Mast Assembly
107918	SA180 Kneepad Outer Pivot Assembly
107920	SA180 Kneepad Inner Pivot Assembly
107308	SA180 Kneepad Assembly
107260	SA180 Foot Tray Assembly
107901	Front Castor Kit (Pair)
107902	Rear Castor Kit (Pair)
510005	Battery Pack (2.0Ah Lead Acid)
510091	Mains Snap Lead
510092	Charger Lead (UK Plug)
107088	Lift Actuator Lead (500mm)
999065	Mackworth SA180 User Manual





PRODUCT VARIATIONS

- 107750 - Mackworth SA180





3. Part Code Detail

	<p>Leg Rod Assembly 107908</p>
	<p>Leg Actuator Lead (750mm) 119105</p>
	<p>4 Way Handset (Single Battery LED) 123105</p>
	<p>4 Way Control Box (Standard) 123108</p>





Part Code Detail

	<p>LH Leg Assembly 107904</p>
	<p>RH Leg Assembly 107903</p>
	<p>Mast Fixing Screw (Small) 000213</p>
	<p>Mast Fixing Screw (Large) 000746</p>





Part Code Detail

 <p>A technical line drawing of the SA180 Base Assembly, showing a rectangular metal base with various mounting points, a vertical support structure, and a circular component on the right side.</p>	<p>SA180 Base Assembly 107212</p>
 <p>A technical line drawing of the Lift Actuator Kit, featuring a long, cylindrical actuator with a handle on the left and a threaded end on the right.</p>	<p>Lift Actuator Kit 107915</p>
 <p>A technical line drawing of the Leg Actuator Kit, showing a rectangular motor housing with a cylindrical shaft extending from the left and a circular component on the right. The text 'DO NOT STAND' is visible on the top of the motor housing.</p>	<p>Leg Actuator Kit 107928</p>
 <p>A technical line drawing of the SA180 Boom Assembly, depicting a long, curved boom with various joints and a circular component at the end.</p>	<p>SA180 Boom Assembly 107921</p>





Part Code Detail

 A technical line drawing of a long, horizontal metal mast assembly. It features several mounting points, including a hook-like end on the right and various brackets and bolts along its length.	<p>SA180 Mast Assembly 107922</p>
 A technical line drawing of a rectangular metal component with a curved bracket on one side. It has a cylindrical end on the right and a square-shaped cutout on the left.	<p>SA180 Kneepad Outer Pivot Assembly 107918</p>
 A technical line drawing of a long, thin metal bar with a cylindrical end on the left and a series of four small rectangular slots along its top edge.	<p>SA180 Kneepad Inner Pivot Assembly 107920</p>
 A technical line drawing of a complex, multi-part assembly consisting of several rectangular metal components joined together, forming a large, rectangular structure with rounded corners and a handle-like protrusion on the right.	<p>SA180 Kneepad Assembly 107308</p>

Part Code Detail

	<p>Foot Tray Assembly 107260</p>
	<p>Front Castor Kit (Pair) 107901</p>
	<p>Rear Castor Kit (Pair) 107902</p>
	<p>Battery Pack (2.0Ah Lead Acid) 510005</p>

Part Code Detail

	<p>Mains Snap Lead 510091</p>
	<p>Charger Lead (UK Plug) 510092</p>
	<p>Lift Actuator Lead (500mm) 107088</p>
<div data-bbox="259 1141 400 1369" style="border: 1px solid black; padding: 5px; text-align: center;">  USER MANUAL </div>	<p>Mackworth SA180 User Manual 999065</p>



4. Notes



Notes

5. Contacts

DEALER CONTACT DETAILS

MANUFACTURER CONTACT DETAILS

Unit 1, Tir Llwyd Industrial Estate, St Asaph Avenue, Kinmel Bay,
Conwy, United Kingdom, LL18 5JZ

Telephone: +44 (0)1924 840100
www.prismmedical.co.uk

DISCLAIMER

While every effort has been made to ensure the accuracy of information contained in this manual, no liability can be accepted by Prism Medical UK for any errors or omissions. Prism Medical UK operates a policy of continuous improvement. Specifications and other data are subject to change without notice.

prism**medicaluk**
a prism**healthcare** company

UK
CA



CE